

MARKERS AND REFRACTIVE

Optical Zone Markers	3 – 1
Corneal Markers.....	3 – 1
Centration Marker	3 – 1
Corneal Transplant Markers	3 – 2
Glaucoma Markers.....	3 – 2
LRI/Toric Markers	3 – 2
LASIK Markers.....	3 – 4
PRK, LASIK, LASEK, IntraLasik and Flap Management Instruments.....	3 – 4
LASIK Fixation Ring	3 – 8
LASIK Cannulas	3 – 8
LASEK Instruments.....	3 – 9
Intralasik Instruments	3 – 11
PRL & ICL Implantation	3 – 12

MARKERS and REFRACTIVE

Bores Axis Marker



AS3-001

- To be used with Mendez Degree Gauge or Grooved Mendez Degree Gauge

OPTICAL ZONE MARKERS

Hoffer Optical Zone Marker With Cross Hairs

AS3-100T-3.0

- 3mm, Titanium

AS3-100T-3.25

- 3.25mm, Titanium

AS3-100T-3.5

- 3.5mm, Titanium

AS3-100T-3.75

- 3.75mm, Titanium

AS3-100T-4.0

- 4mm, Titanium

AS3-100T-4.25

- 4.25mm, Titanium

AS3-100T-4.5

- 4.5mm, Titanium

AS3-100T-4.75

- 4.75mm, Titanium

AS3-100T-5.0

- 5mm, Titanium

AS3-100T-5.5

- 5.5mm, Titanium

AS3-100T-6.0

- 6mm, Titanium

AS3-100T-7.0

- 7mm, Titanium

AS3-100T-8.0

- 8mm, Titanium

AS3-100T-9.0

- 9mm, Titanium



CORNEAL MARKERS

Osher-Neuman Low Profile Corneal Marker



AS3-200T-4

- 4 Lines, Titanium

AS3-200T-6

- 6 Lines, Titanium

AS3-200T-8

- 8 Lines, Titanium

AS3-200T-10

- 10 Lines, Titanium

AS3-200T-12

- 12 Lines, Titanium

AS3-200T-16

- 16 Lines, Titanium

CENTRATION MARKERS

Whipple Capsulorhexis Centration Marker



AS3-300S

- Marks 5.5mm-6mm
- Stainless Steel

MARKERS and REFRACTIVE

CORNEAL TRANSPLANT MARKERS

Corneal Transplant Marker



AS3-400

- Internal Diameter 7mm
- External Diameter 11mm



AS3-401S

DSEK MARKERS

Corneal DSEK



AS3-410-8.0

- 8.0mm

AS3-410-9.0

- 9.0mm

AS3-410-8.5

- 8.5mm

AS3-410-9.5

- 9.5mm

Corneal DSEK Double Ended



AS3-411

- 8.5mm
- 9.5mm Marker with Cross Hairs

LRI/TORIC MARKERS

LRI Marker, 6mm & 8mm Chord Lengths



AS3-500

- 6mm & 8mm Chord Lengths With Degree Gauge (Uses Dell/Gills Nomograms)

LRI Marker, 40-60-80 Degrees



AS3-501

- 40-60-80 degrees (Uses Nichamin Nomograms)

MARKERS and REFRACTIVE

Toric IOL Marker



AS3-502

LRI Slit Lamp Marker



AS3-503

- To Be Used With Mendez Gauge
- Used To Mark Horizontal Axis

LRI Slit Lamp Marker - 3 Points



AS3-504

- Used To Mark Horizontal Axis
- 3 Points

LRI Slit Lamp Gravity Marker



AS3-505

- Used To Mark Horizontal Axis

LRI/Toric Gravity Marker



AS3-506

- 2 Points
- Used To Mark Horizontal Axis

LRI/Toric IOL Slit Lamp Gravity Marker



AS3-507

- 3 Points

MARKERS and REFRACTIVE

Mendez Degree Gauge With Fixation




AS1-355

- Atraumatic Teeth For Fixation
- Calibrated Marks Every 10 Degrees From 0-180
- Internal 13mm / External 19mm

LASIK MARKERS

LASIK Marker



AS3-600

- Chayet Style
- For Identification and Realignment of Lamellar Flap

Lasik Marker, with 2 Radials and 1 Para Radials



AS3-601

- Lu Style
- For Identification and Realignment of Lamellar Flap

Lasik Marker, with 3.5mm and 4.0mm Intersecting Circles



AS3-602

- Titanium (Gulani Style)
- For Identification and Realignment of Lamellar Flap

PRK, LASIK, LASEK, INTRALASIK AND FLAP MANAGEMENT INSTRUMENTS

Paton Double Ended Spatula And Spoon



AS8-130

- For Corneal Button Transfer

MARKERS and REFRACTIVE

Hockey Epithelium Removal Knife



AS4-022

- Sharp 2.3mm Front Edge, Semi-Sharp 4mm Bottom Edge

AS4-023

- Sharp 2.7mm Front Edge, Semi-Sharp 6mm Bottom Edge

Lindstrom LASIK/PRK Spatula And Epithelium Removal Board



AS3-607

- Semi-Sharp end of Instrument removes epithelium
- Fine Spatula used for Flap Manipulation during LASIK Surgery

Bahamond Retreatment LASIK Spatula



AS3-608

Bloomberg Epithelial Debrider



AS3-609

- 5mm x 2mm Semi-Sharp Spatulated Tip

LASIK/PRK Spatula and Epithelium Removal



AS3-611

- Lindstrom Style

MARKERS and REFRACTIVE

Buratto Flap Protector



AS3-612

- Nasal Hinge

LASIK Flap Smoother



AS3-613

LASIK Flap Manipulator

AS3-620

Machat PRK/LASIK Flap Spatula

AS3-621

- 2mm x 5mm Golf Club Shaped Blade with Semi-Sharp Outer Edges

LASIK Flap Spatula



AS3-622

- Double Ended
- Banaji Style

MARKERS and REFRACTIVE

LASIK Flap Conformer



AS3-623

LASIK Spatula And Flap Retreatment Instrument



AS3-624

- Bahamond

LASIK Retreatment Spatula



AS3-625

- Bahamond

LASIK Retreatment Separator

AS3-626

LASIK Flap Separator

AS3-627

- Katzin
- Double Ended



AS3-628

MARKERS and REFRACTIVE

Disposable Marking Pen



AS25-001
• Pack of 10

Disposable Marking Pad



AS25-002
• Gentian Violet Pack of 10

LASIK FIXATION RING

LASIK Hinge Protector Fixation Ring



AS1-353
• 11mm Inner Diameter
• Swivel
• Titanium

LASIK CANNULAS

Vidaurri LASIK Double Cannula



AS13-1700-H
• 25 Gauge
(Vidaurri-Horizontal)

AS13-1700-V
• 25 Gauge
(Vidaurri-Vertical)

LASIK Irrigation Cannula



AS13-1701-25g
• 4 Ports
• Curved
• 25 Gauge

AS13-1701-27g
• 4 Ports
• Curved
• 27 Gauge

MARKERS and REFRACTIVE

LASIK Cannula, Slade Type

**AS13-1703**

- Spatulated Tip Opening
- 26 Gauge

LASIK Irrigation Cannula, Buratto Style

**AS13-1704-23g**

- 4 Ports
- Curved
- 25 Gauge

AS13-1704-27g

- 4 Ports
- Curved
- 27 Gauge

LASIK Fountain Cannula, Gimbel Style

**AS13-1705-L**

- Left
- 25 Gauge

AS13-1705-R

- Right
- 25 Gauge

LASEK INSTRUMENTS

LASEK Trephine

**AS3-700**

- 8.0mm
- Preset Depth
70 Microns

AS3-701

- 9.0mm
- Preset Depth
70 Microns

LASEK Funnel

**AS3-702**

- 8.5 mm

AS3-703

- 9.5 mm

MARKERS and REFRACTIVE

LASEK 8mm Trephine



AS3-704

- 8.5mm Funnel

AS3-705

- 9.5mm Funnel

LASEK Detaching Spatula



AS3-706

LASEK Epithelial Micro Hoe



AS3-707

- Bahamon

LASEK Epithelial Flap Repositioning Spatula



AS3-708

LASEK Bow Dissector



AS3-709

MARKERS and REFRACTIVE

LASEK Spatula



AS3-710

LASEK Knife



AS3-711

LASEK Alcohol Cannula



AS13-1720



AS13-1721

• Disk Shape

INTRALASIK INSTRUMENTS

Intralasik Flap Spatula Double Ended



AS3-800

MARKERS and REFRACTIVE

Intralask Flap Spatula



AS3-801

PRL & ICL IMPLANTATION

ICL Cartridge Loading Forceps



AS3-900
• Stainless Steel

PRL Implantation Forceps

AS3-901
• Angled
• Titanium

Forceps for Artisan Phakic Lens



AS3-902

ICL Removal Forceps



AS3-904
• Titanium

MARKERS and REFRACTIVE

Coaxial Angled ICL Loading Forceps

**AS3-905**

• Tip Only

Pallikaris ICL Manipulator

**AS3-930**

Battle ICL Manipulator

**AS3-931**

Deitz ICL Tucker and Slider

**AS3-932**

