

SAPPHIRE AND DIAMOND KNIVES

SAPPHIRE KNIVES

Sapphire Universal & Side Port Knives	6 – 1
Sapphire Crescent Tunnel Knives	6 – 2
Sapphire Clear Cornea Phaco Spear Knives	6 – 2
Sapphire Trapezoid Knives	6 – 2

BLACK DIAMOND KNIVES

Black Diamond Universal & Side Port Knives....	6 – 3
Black Diamond Crescent Tunnel Knives	6 – 3
Black Diamond Clear Cornea Phaco Slit Knives	6 – 3

DIAMOND KNIVES

Universal Straight Knives	6 – 4
Angled Phaco Knives	6 – 4
Crescent Knives	6 – 5
Trapezoid Self-Diving Knives.....	6 – 5
Side Port & MICS Knives	6 – 5
Step Knives for Cataract Surgery and LRI.....	6 – 6

SAPPHIRE KNIVES

SAPPHIRE UNIVERSAL & SIDE PORT KNIVES

Sapphire Double Edge



AS6-001

- 1.0mm
- 30 degrees
- Straight Handle



AS6-002

- 1.0mm
- 45 degrees
- Straight Handle



Sapphire Lancet Blade

AS6-003

- .50mm Blade
- Straight Handle



AS6-004

- 1.0mm Blade
- Straight Handle

Sapphire Straight Double Lancet

AS6-005

- 60 degree Angle
- .75mm Blade
- Straight Handle



AS6-006

- 60 degree Angle
- 1.0mm Blade
- Straight Handle

Sapphire Trifacet Blade

AS6-007

- 1.0mm
- Straight Handle



AS6-008

- 1.2mm
- Straight Handle

Sapphire Zaldiver Blade


AS6-009

- Straight Handle



SAPPHIRE KNIVES

SAPPHIRE CRESCENT TUNNEL KNIVES

**AS6-020**

- 1.5mm wide
- Angle Handle

AS6-022

- 2.2mm wide
- Angled Handle

AS6-021

- 2.0mm wide
- Angled Handle

AS6-023

- 2.5mm wide
- Angled Handle

SAPPHIRE CLEAR CORNEA PHACO SPEAR KNIVES

**AS6-030**

- 1.0mm
- 60 degrees
- Angled Handle

AS6-031

- 1.2mm
- 60 degrees
- Angled Handle

AS6-032

- 1.5mm
- 60 degrees
- Angled Handle

AS6-033

- 2.2mm
- 60 degrees
- Angled Handle

AS6-034

- 2.5mm
- 60 degrees
- Angled Handle

AS6-035

- 2.65mm
- 60 degrees
- Angled Handle

AS6-036

- 2.8mm
- 60 degrees
- Angled Handle

AS6-037

- 3.0mm
- 60 degrees
- Angled Handle

AS6-038

- 3.2mm
- 60 degrees
- Angled Handle

SAPPHIRE TRAPEZOID KNIVES

**AS6-040**

- 1.8-2.2mm Blade
(All Sharp)

AS6-041

- 2.0-2.3mm Blade
(All Sharp)

AS6-042

- 2.2-2.5mm Blade
(All Sharp)

AS6-043

- 2.5-2.9mm Blade
(All Sharp)

AS6-044

- 2.55-3.0mm
Blade (All Sharp)

BLACK DIAMOND KNIVES

BLACK DIAMOND UNIVERSAL & SIDE PORT KNIVES

Black Diamond Double Edge



- AS6-075**
- 1.0mm
 - 30 degrees
 - Straight Handle



- AS6-076**
- 1.0mm
 - 45 degrees
 - Straight Handle

Black Diamond Straight Double Lancet



- AS6-077**
- 35 degree Angle
 - 1.0mm Blade
 - Straight Handle

- AS6-078**
- 60 degree Angle
 - 1.0mm Blade
 - Straight Handle

- AS6-079**
- 60 degree Angle
 - 1.2mm Blade
 - Straight Handle

BLACK DIAMOND CRESCENT TUNNEL KNIVES




- AS6-080**
- 1.5mm wide
 - Angled Handle

- AS6-081**
- 2.0mm wide
 - Angled Handle

- AS6-082**
- 2.5mm wide
 - Angled Handle

- AS6-083**
- 2.8mm wide
 - Angled Handle



BLACK DIAMOND CLEAR CORNEA PHACO SLIT KNIVES



- AS6-091**
- 2.5mm
 - 60 degrees
 - Angled Handle

- AS6-092**
- 2.65mm
 - 60 degrees
 - Angled Handle

- AS6-093**
- 2.8mm
 - 60 degrees
 - Angled Handle

- AS6-094**
- 3.0mm
 - 60 degrees
 - Angled Handle

- AS6-095**
- 3.2mm
 - 60 degrees
 - Angled Handle

DIAMOND KNIVES

UNIVERSAL STRAIGHT KNIVES

30 Degree Single Edge Blade



AS6-100
 • Straight Titanium
 Universal Handle

45 Degree Single Edge Blade



AS6-101
 • Straight Titanium
 Universal Handle



45 Degree Double Edge Blade



AS6-102
 • Straight Titanium
 Universal Handle

Trifacet Blade



AS6-103
 • 1.0mm
 • Straight Titanium
 Universal Handle

20 Degree Trifacet Blade



AS6-104
 • .2mm-1.0mm Width
 • Straight Titanium
 Universal Handle

Lancet Blade



AS6-105
 • 1.0mm
 • Straight Titanium
 Universal Handle

Brown-Osher Five Facet Diamond Blade



AS6-106
 • Straight Titanium
 Universal Handle

ANGLED PHACO KNIVES

Clear Cornea Blades



AS6-120
 • 2.2mm
 • Angled Titanium
 Universal Handle

AS6-121
 • 2.5mm
 • Angled Titanium
 Universal Handle

AS6-122
 • 2.8mm
 • Angled Titanium
 Universal Handle

AS6-123
 • 3.0mm
 • Angled Titanium
 Universal Handle

AS6-124
 • 3.2mm
 • Angled Titanium
 Universal Handle

DIAMOND KNIVES

CRESCENT KNIVES

**AS6-130**

- 1.5mm Round
- Angled Titanium Universal Handle

AS6-131

- 2.0mm Round
- Angled Titanium Universal Handle

TRAPEZOID SELF-DIVING DIAMOND KNIVES

**AS6-140**

- 1.8 / 2.0mm
- Angled Titanium Handle

AS6-141

- 2.0 / 2.3mm
- Angled Titanium Universal Handle

AS6-142

- 1.9 / 2.5mm
- Angled Titanium Universal Handle

AS6-143

- 2.3 / 2.8mm
- Angled Titanium Universal Handle

AS6-144

- 2.6 / 3.0mm
- Angled Titanium Universal Handle

AS6-145

- 2.7 / 3.2mm
- Angled Titanium Universal Handle

SIDE PORT & MICS KNIVES

Universal ICL Knife

**AS6-150**

- .55mm-1.0mm Width
- Angled Titanium Universal Handle

Trapezoid Side Port Knife

**AS6-151**

- 1.0mm Width
- Angled Titanium Handle

Micro Incision Knife

**AS6-152**

- 0.8 / 1.2 mm Width
- Angled Titanium Universal Handle

AS6-153

- 1.3 / 1.5mm Width
- Angled Titanium Universal Handle

AS6-154

- 1.5 / 1.8mm Width
- Angled Titanium Universal Handle

AS6-155

- 1.8 / 2.0mm Width
- Angled Titanium Universal Handle

AS6-156

- 0.9 / 1.3mm Width
- Angled Titanium Universal Handle

STEP KNIVES FOR CATARACT SURGERY AND LRI

Multi-Step Knife



AS6-200

- 1.0mm Width
- 0.50, 0.55 & 0.60mm, & Full Extension
- Step Handle with Footplate



Multi-Incision 5 Facet Blade



AS6-201

- 1.0mm Width
- 0.50, 0.55 & 0.60mm, & Full Extension
- Step Handle with Footplate

Single-Step Knives



AS6-202

- 1.0mm Width
- Preset Depth 0.50mm

AS6-203

- 1.0mm Width
- Preset Depth 0.60mm



AS6-202-SL

- 1.0mm Width
- Preset Depth 0.50mm
- LRI Slit Lamp

AS6-203-SL

- 1.0mm Width
- Preset Depth 0.60mm
- LRI Slit Lamp

Diamond Knife Cleaning Pack



AS25-005

- 1-Pack
- The pack consists of 3 solutions that are used to clean residual debris from blades or instruments prior to sterilization.